

For Research Use Only

# PE Anti-Mouse CD115 (c-fms) (AFS98)

Catalog Number: PE-65095



## Basic Information

Catalog Number:

PE-65095

Size:

100ug , 0.2 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC036343

GeneID (NCBI):

12978

UNIPROT ID:

P09581

Full Name:

colony stimulating factor 1 receptor

Purification Method:

Affinity purification

CloneNo.:

AFS98

Excitation/Emission maxima  
wavelengths:

496 nm, 565 nm / 578 nm

## Applications

Tested Applications:

FC

Species Specificity:

mouse

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> Balb/c mouse peritoneal macrophages were surface stained with 0.125 ug PE Anti-Mouse CD115 (c-fms) (AFS98) (PE-65095, Clone: AFS98) (red), or 0.125 ug PE Rat IgG2a Isotype Control (2A3) (PE-65209, Clone: 2A3) (blue). Cells were not fixed.