

For Research Use Only

PE Anti-Mouse LPAM-1 (DATK32)

Catalog Number: PE-65101



Basic Information

Catalog Number: PE-65101	GenBank Accession Number: BC011184	Purification Method: Affinity purification
Size: 100ug, 0.2 mg/ml	GeneID (NCBI): 16421	CloneNo.: DATK32
Source: Rat	Full Name: integrin beta 7	Excitation/Emission maxima wavelengths: 496 nm, 565 nm / 578 nm
Isotype: IgG2a, kappa		

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 mouse bone marrow cells were surface co-stained with CoraLite® Plus 488 Anti-Mouse CD45R (B220) (RA3-6B2) and 0.2 μ g PE Anti-Mouse LPAM-1 (PE-65101, Clone:DATK32) or 0.2 μ g Isotype Control. Cells were not fixed.