

For Research Use Only

# PE Anti-Mouse LPAM-1 (DATK32)

Catalog Number: PE-65101



## Basic Information

Catalog Number:

PE-65101

Size:

100ug, 0.2 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC011184

GeneID (NCBI):

16421

UNIPROT ID:

P26011

Full Name:

integrin beta 7

Purification Method:

Affinity purification

CloneNo.:

DATK32

Excitation/Emission maxima wavelengths:

496 nm, 565 nm / 578 nm

## Applications

Tested Applications:

FC

Species Specificity:

mouse

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> mouse bone marrow cells were surface co-stained with CoraLite® Plus 488 Anti-Mouse CD45R (B220) (RA3-6B2) and 0.2 ug PE Anti-Mouse LPAM-1 (PE-65101, Clone:DATK32) or 0.2 ug Isotype Control. Cells were not fixed.