

For Research Use Only

# PE Anti-Mouse CD69 (H1.2F3)

Catalog Number: PE-65105

1 Publications



## Basic Information

<b>Catalog Number:</b> PE-65105	<b>GenBank Accession Number:</b> BC106997	<b>Purification Method:</b> Affinity purification
<b>Size:</b> 100ug, 0.2 mg/ml	<b>GeneID (NCBI):</b> 12515	<b>CloneNo.:</b> H1.2F3
<b>Source:</b> Armenian Hamster	<b>Full Name:</b> CD69 antigen	<b>Excitation/Emission maxima wavelengths:</b> 496 nm, 565 nm / 578 nm
<b>Isotype:</b> IgG		

## Applications

**Tested Applications:**  
FC

**Cited Applications:**  
FC

**Species Specificity:**  
Mouse

**Cited Species:**  
mouse

## Background Information

CD69, also known as AIM, EA-1, Leu-23, and MLR3, is a type II transmembrane glycoprotein that belongs to the C-type lectin superfamily (PMID: 8340758; 7804122). CD69 is constitutively expressed by mature thymocytes, platelets, several subsets of tissue resident immune cells (including resident memory T cells and gamma delta T cells), and is inducibly expressed by activated T cells, B cells, natural killer (NK) cells, monocytes, neutrophils (PMID: 8100776; 28475283). CD69 has been identified as an early activation marker of lymphocytes and is commonly used as a marker of activated lymphocytes and NK cells (PMID: 28475283; 25759842). It is involved in the regulation of immune responses (PMID: 15745855).

## Notable Publications

Author	Pubmed ID	Journal	Application
Linxiang Cai	37578941	J Med Chem	FC

## Storage


**Storage:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.1% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> unstimulated or anti-CD3/CD28 treated mouse splenocytes were surface stained with 0.25 ug PE Anti-Mouse CD69 (PE-65105, Clone:H1.2F3). Cells were not fixed.