

For Research Use Only

# PE Anti-Mouse CD24 (M1/69)

Catalog Number: PE-65121

1 Publications



## Basic Information

Catalog Number:

PE-65121

Size:

100ug, 0.2 mg/ml

Source:

Rat

Isotype:

IgG2b, kappa

GenBank Accession Number:

BC075622

GeneID (NCBI):

12484

UNIPROT ID:

P24807

Full Name:

CD24a antigen

Purification Method:

Affinity purification

CloneNo.:

M1/69

Excitation/Emission maxima  
wavelengths:

496 nm, 565 nm / 578 nm

## Applications

Tested Applications:

FC

Cited Applications:

FC

Species Specificity:

Mouse

Cited Species:

mouse

## Notable Publications

Author	Pubmed ID	Journal	Application
Yao Qi	39603539	J Control Release	FC

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> BALB/c mouse splenocytes were surface stained with 0.06 ug PE Anti-Mouse CD24 (PE-65121, Clone: M1/69) (red) or 0.06 ug isotype control antibody (black). Cells were not fixed.