

For Research Use Only

PE Anti-Mouse CD24 (M1/69)

Catalog Number: PE-65121



Basic Information

Catalog Number: PE-65121	GenBank Accession Number: BC075622	Purification Method: Affinity purification
Size: 100ug, 0.2 mg/ml	GeneID (NCBI): 12484	CloneNo.: M1/69
Source: Rat	Full Name: CD24a antigen	Excitation/Emission maxima wavelengths: 496 nm, 565 nm / 578 nm
Isotype: IgG2b, kappa		

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ BALB/c mouse splenocytes were surface stained with 0.06 ug PE Anti-Mouse CD24 (PE-65121, Clone: M1/69) (red) or 0.06 ug isotype control antibody (black). Cells were not fixed.