

For Research Use Only

PE Anti-Mouse CD197 (4B12)

Catalog Number: PE-65158



Basic Information

Catalog Number:
PE-65158

Size:
100ug, 0.5 mg/ml

Source:
Rat

Isotype:
IgG2a kappa

GenBank Accession Number:

GeneID (NCBI):
12775

Full Name:
chemokine (C-C motif) receptor 7

Purification Method:

The purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. The conjugate is purified by size-exclusion chromatography.

CloneNo.:
4B12

Excitation/Emission maxima
wavelengths:
496 nm, 565 nm / 578 nm

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

Storage

Storage:
Store at 2-8°C. Avoid exposure to light.
Storage Buffer:
PBS with 0.1% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ C57BL/6 mouse splenocytes were surface stained with 2.50 ug PE Anti-Mouse CD197 (PE-65158, Clone: 4B12) (red) or 2.50 ug isotype control antibody (blue). Cells were not fixed.