

For Research Use Only

# PE Anti-Human CD120a (H398)

Catalog Number: PE-65170



## Basic Information

Catalog Number:

PE-65170

Size:

100 tests, 10 µl/test

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

BC010140

GeneID (NCBI):

7132

Full Name:

tumor necrosis factor receptor  
superfamily, member 1A

Calculated MW:

455 aa, 50 kDa

Purification Method:

The purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. The conjugate is purified by size-exclusion chromatography.

CloneNo.:

H398

Excitation/Emission maxima  
wavelengths:

496 nm, 565 nm / 578 nm

## Applications

Tested Applications:

FC

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



100 ul human peripheral blood were surface stained with 10.00 ul PE Anti-Human CD120a (PE-65170, Clone: H398). Cells were not fixed.



100 ul human peripheral blood were surface stained with PE-mouse IgG2a isotype control antibody. Cells were not fixed.