

For Research Use Only

# PE Anti-Human CD41a (HIP8)

Catalog Number: PE-65173



## Basic Information

Catalog Number:

PE-65173

Size:

100 tests, 5 ul/test

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

BC117443

GeneID (NCBI):

3674

UNIPROT ID:

P08514

Full Name:

integrin, alpha 2b (platelet glycoprotein IIb of IIb/IIIa complex, antigen CD41)

Purification Method:

Unknown

CloneNo.:

HIP8

Excitation/Emission maxima wavelengths:

496 nm, 565 nm / 578 nm

## Applications

Tested Applications:

FC

Species Specificity:

human, non-human primates

## Background Information

The integrins are a superfamily of cell adhesion receptors that bind to extracellular matrix ligands, cell-surface ligands, and soluble ligands (PMID: 17543136). They are transmembrane  $\alpha\beta$  heterodimers and at least 24 distinct integrin heterodimers are formed by the combination of 18  $\alpha$  and eight  $\beta$  known subunits (PMID: 17543136; 20029421). In addition to mediating cell adhesion, integrins also play important roles in modulating signal transduction pathways that control cellular responses including migration, proliferation, differentiation, and apoptosis (PMID: 19118207). Integrin alpha-IIb (ITGA2B, CD41) is expressed on platelets, megakaryocytes and some hematopoietic progenitor cells (PMID: 11934866).

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



100 ul human platelet-rich plasma (diluted 1:2 with PBS) was surface stained with 5 ul PE Anti-Human CD41a (PE-65173, Clone:HIP8) or PE Mouse IgG1 Isotype Control (PE-65124). Cells were not fixed.



100 ul human platelet-rich plasma (diluted 1:2 with PBS) was surface stained with 5 ul PE Anti-Human CD41a (PE-65173, Clone:HIP8) or PE Mouse IgG1 Isotype Control (PE-65124). Cells were not fixed.