

For Research Use Only

PE Anti-Human IgD (IA6-2)

Catalog Number: PE-65175



Basic Information

Catalog Number:
PE-65175

Size:
100tests, 10 µl/test

Source:
Mouse

Isotype:
IgG2a, kappa

GenBank Accession Number:
BC021276

GeneID (NCBI):
3495

Full Name:
immunoglobulin heavy constant
delta

Calculated MW:
573 aa, 63 kDa

Purification Method:

The purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. The conjugate is purified by size-exclusion chromatography.

CloneNo.:
IA6-2

Excitation/Emission maxima wavelengths:
496 nm, 565 nm / 578 nm

Applications

Tested Applications:
FC

Species Specificity:
Human

Background Information

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.1% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human peripheral blood lymphocytes were surface stained with APC Anti-Human CD19 (APC-65110, Clone:HI B19) and 10.00 ul/test PE Anti-Human IgD (PE-65175, Clone:IA6-2). Sample were not fixed.



1X10⁶ human peripheral blood lymphocytes were surface stained with APC Anti-Human CD19 (APC-65110, Clone:HI B19) and PE mouse IgG2a isotype control antibody. Sample were not fixed.