

For Research Use Only

PE Mouse IgG2a Isotype Control (MOPC-173)



Catalog Number: PE-65244

Basic Information

Catalog Number:
PE-65244

Size:
100ug, 0.2 mg/mL

Source:
Mouse

Isotype:
IgG2a, kappa

GenBank Accession Number:
GeneID (NCBI):

Full Name:

Purification Method:
Affinity purification

CloneNo.:
MOPC-173

Excitation/Emission maxima
wavelengths:
496 nm, 565 nm / 578 nm

Applications

Tested Applications:
FC

Species Specificity:

Background Information

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human PBMCs were surface stained with 0.2 ug PE Anti-Mouse IgG2a Isotype Control (PE-65244, Clone:MOPC-173) (red), PE Anti-Human CD45RO (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.