

For Research Use Only

# PE Anti-Mouse CD45 Rabbit Recombinant Antibody

Catalog Number: PE-98035



## Basic Information

Catalog Number:

PE-98035

Size:

100ug , 200 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001111316.2

GeneID (NCBI):

19264

UNIPROT ID:

P06800-4

Full Name:

protein tyrosine phosphatase,  
receptor type, C

Calculated MW:

145 kDa

Purification Method:

Protein A purification

CloneNo.:

240356D1

Excitation/Emission maxima  
wavelengths:

496 nm, 565 nm / 578 nm

## Applications

Tested Applications:

FC

Species Specificity:

mouse

## Background Information

CD45, also known as protein tyrosine phosphatase, receptor type C, is a type I transmembrane protein expressed on the surface of all haematopoietic cells with the exception of erythrocytes and platelets (PMID: 3489673; 28615666). CD45 is a pan-haematopoietic cell marker and has been shown to be essential for T- and B-cell activation and signalling (PMID: 9429890; 16378097).

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> mouse splenocytes were surface stained with 0.25 ug PE Anti-Mouse CD45 Rabbit Recombinant Antibody (PE-98035, Clone: 240356D1) (red) or Isotype Control (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were not fixed.



1x10<sup>6</sup> mouse splenocytes were surface stained with 0.25 ug PE Anti-Mouse CD45 Rabbit Recombinant Antibody (PE-98035, Clone: 240356D1) or Isotype Control, and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were not fixed.