

For Research Use Only

# FcZero-rAb® PE Anti-Human CD28 Rabbit Recombinant Antibody

Catalog Number: PE-FcA98032



## Basic Information

**Catalog Number:**

PE-FcA98032

**Size:**

100 tests, 5 µl/test

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

EG0889

**GenBank Accession Number:**

BC093698

**GeneID (NCBI):**

940

**ENSEMBL Gene ID:**

ENSG00000178562

**UNIPROT ID:**

P10747

**Full Name:**

CD28 molecule

**Calculated MW:**

220 aa, 25 kDa

**Purification Method:**

Protein A purification

**CloneNo.:**

240354D2

**Recommended Dilutions:**

FC: 5 µl per 10<sup>6</sup> cells in a 100 µl suspension

**Excitation/Emission maximum wavelengths:**

496 nm, 565 nm / 578 nm

## Applications

**Tested Applications:**

FC

**Species Specificity:**

human, monkey, non-human primates

**Positive Controls:**

FC : human PBMCs, cynomolgus PBMCs

## Background Information

CD28 (T-cell-specific surface glycoprotein CD28), also known as T44 and Tp44, is a 44 kD disulfide-linked homodimeric type I glycoprotein (PMID: 2162180). It is a member of the immunoglobulin superfamily and is expressed on most T lineage cells, NK cell subsets, and plasma cells (PMID: 2162180, 8386518). CD28 may affect in vivo immune responses by functioning both as a cell adhesion molecule linking B and T lymphocytes and as the surface component of a novel signal transduction pathway (PMID: 2162180, 3021470). CD28 binds both CD80 and CD86 with a highly conserved motif MYPPY in the CDR3-like loop (PMID: 15696168, 7964482). CD28 is considered a major co-stimulatory molecule, inducing T lymphocyte activation and IL-2 synthesis, and preventing cell death (PMID: 1348520).

## Storage

**Storage:**

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data




1x10<sup>6</sup> human PBMCs were surface stained with CL750-FcA65570, and 5 ul PE Anti-Human CD28 Rabbit RecAb (PE-FcA98032, Clone: 240354D2) or PE Rabbit IgG Isotype Control RecAb (PE-FcA98136, Clone: 240953C9). Cells were not fixed. Lymphocytes were gated.



1x10<sup>6</sup> human PBMCs were surface stained with 5 ul PE Anti-Human CD28 Rabbit RecAb (PE-FcA98032, Clone: 240354D2) (red) or PE Rabbit IgG Isotype Control RecAb (PE-FcA98136, Clone: 240953C9) (blue). Cells were not fixed. Lymphocytes were gated.



1x10<sup>6</sup> cynomolgus PBMCs were surface stained with 5 ul PE Anti-Human CD28 Rabbit RecAb (PE-FcA98032, Clone: 240354D2) or PE Rabbit IgG Isotype Control RecAb (PE-FcA98136, Clone: 240953C9). Cells were not fixed.



1x10<sup>6</sup> cynomolgus PBMCs were surface stained with 5 ul PE Anti-Human CD28 Rabbit RecAb (PE-FcA98032, Clone: 240354D2) (red) or PE Rabbit IgG Isotype Control RecAb (PE-FcA98136, Clone: 240953C9) (blue). Cells were not fixed.