

For Research Use Only


Foxp3 / Transcription Factor Fix/Perm Concentrate (4X)

Catalog Number: PF00011-A



Description	Foxp3 / Transcription Factor Fix/Perm Concentrate (4X): provided in a concentrated soltion which, when diluted with Foxp3 / Transcription Factor Fix/Perm Diluent (1X) provide best results during intracellular staining of transcription factors using fluorescently conjugated antibodies.
Package	50 mL
Storage	4°C
Expiration Date	12 months after receiving.

Validation Data



1X10⁶ mouse splenocytes were surface stained with CoraLite488-conjugated Anti-Mouse CD4 (GK1.5) (CL488-65104, Clone: GK1.5) and then fixed with 1X Transcription Factor Fix/Perm buffer (PF00011-A) and permeabilized with 1X Flow Cytometry Perm Buffer (PF00011-C). Cells were then stained with CoraLite647-conjugated mouse... IgG1 isotype control or 5 ul CoraLite647-conjugated Anti-Mouse Foxp3 (CL647-65089, Clone: 3G3).

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB(1-888-478-4522)(toll free in USA), or
1(312) 455-8498(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under
Proteintech Group brand and is not available
to purchase from any other manufacturer.