

For Research Use Only.

Foxp3 / Transcription Factor Fix/Perm Diluent (1X)

Catalog Number: PF00011-B



Description

Foxp3 / Transcription Factor Fix/Perm Diluent (1X) Intended for use as a diluent for Foxp3 / Transcription Factor Fix/Perm Concentrate (4X)

Package

160 mL

Storage

4°C

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB(1-888-478-4522)(toll free in USA), or 1(312) 455-8498(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.