

For Research Use Only

# Exosome Isolation and Protein Extraction Kit (Cell Culture Supernatant/Urine)

Catalog Number: PK10028



## Product Information

Component	Size	Storage
Exosome Isolation Buffer (PK10028-A)	20 mL	-20°C
Exosome Lysis Buffer (PK10028-B)	10 mL	-20°C
SDS-PAGE Protein Loading Buffer (4×) (PK10028-C)	2 mL	-20°C
Standard Protease Inhibitor Mixture (100×) (PK10028-D)	0.5 mL	-20°C

## Package

40 T

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB(1-888-478-4522)(toll free in USA), or  
1(312) 455-8498(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.