

For Research Use Only

ZymeBright™ APC Site-Specific Antibody Conjugation Kit

Catalog Number: PK30081



Product Information

ZymeBright™ APC Site-Specific Antibody Conjugation Kit is enzymatic method for the conjugation of APC to the Fc region of Rabbit IgG, Rat IgG1 & IgG2a, Human IgG1-4, or Mouse IgG1 & IgG3 for the use in fluorescent immunoassays. The expected degree of labeling (F/P ratio) of antibodies conjugated with this kit is 2.

Product Components

Cat No.	Name	Size
Component A (PK3007X-A)	Antibody Modifier Reagent (Vial A)	5 uL
Component B (PK30081-B)	Lyophilized Modified APC (Vial B)	1 rxn
Component C (PK3007X-C)	Lyophilized PTGase (Vial C)	3 mg (provided in excess)
Component D (PK3007X-D)	Conjugation Buffer (Vial D)	200 uL
Component E (PK3007X-E)	Quenching Reagent (Vial E)	5 mg
Component F (PK3007X-F)	0.6 mL Low Adhesion Tube	0.6 mL*2
Labeling Scale: 100 ug unconjugated antibody		

Package

100 ug

Storage

Store at -20°C for up to 1 year. After opening, store at 2-8°C. Protect from light and moisture. Reagents must be used within 1 week of opening.


Buffer Compatibility

Buffer Components & Conditions	
Purified antibody	Yes
Antibody in ascites fluid, serum, hybridoma or tissue culture media	No
Antibody concentration	0.5-2 mg/mL
pH	6.5-8.0
Amine free buffer (e.g. MES, MOPS, HEPES, PBS)	Yes
Non-buffering salts (e.g. sodium chloride)	Yes
BSA	Yes
Sodium Azide	Yes
Chelating agents (e.g. EDTA)	Yes
Glycerol	<10%
Sugars	Yes
Gelatin	Unknown
Tris	Yes
Glycine	No
Thiomersal / Thimerosal	Yes
Merthiolate	Yes
Proclin	Yes
Borate buffer	Yes
Primary amines (e.g. lysine, peptide or ethanolamine)	No
Heavy metal ions (Zn ²⁺ , Cu ²⁺ , Mn ²⁺ , or Fe ³⁺)	No

Required Materials

Orbital Thermal Mixer, Pipettes, Sterile PBS.

Validation Data



1x10⁶ human PBMCs were surface stained with 0.1ug APC Anti-Human CD7 Rabbit RecAb (APC-98204, Clone: 242402H8), and 0.1ug FcZero-rAb® CoraLite® Plus 488 Anti-Human CD4 Rabbit Recombinant Antibody (CL488-FcA98042, Clone: 240427E12), and 0.1ug Rabbit IgG Isotype Control RecAb (APC-98136, Clone: 240953C9). Cells were... incubated with MonoZero® Human Fc Blocking Reagent (PF00032) prior to staining. Cells were not fixed.

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB(1-888-478-4522)(toll free in USA), or
1(312) 455-8498(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under
Proteintech Group brand and is not available
to purchase from any other manufacturer.