

For Research Use Only

ZymeBright™ APC-Cyanine7 Site-Specific Antibody Conjugation Kit

Catalog Number: PK30089



Product Information

ZymeBright™ APC-Cyanine7 Site-Specific Antibody Conjugation Kit is enzymatic method for the conjugation of APC-Cyanine7 to the Fc region of Rabbit IgG, Rat IgG1 & IgG2a, Human IgG1-4, or Mouse IgG1 & IgG3 for the use in fluorescent immunoassays. The expected degree of labeling (F/P ratio) of antibodies conjugated with this kit is 2.

Product Components

Cat No.	Name	Size
Component A (PK3007X-A)	Antibody Modifier Reagent (Vial A)	5 uL
Component B (PK30089-B)	Lyophilized Modified APC-Cyanine7 (Vial B)	1 rxn
Component C (PK3007X-C)	Lyophilized PTGase (Vial C)	3 mg (provided in excess)
Component D (PK3007X-D)	Conjugation Buffer (Vial D)	200 uL
Component E (PK3007X-E)	Quenching Reagent (Vial E)	5 mg
Component F (PK3007X-F)	0.6 mL Low Adhesion Tube	0.6 mL*2
Labeling Scale: 100 ug unconjugated antibody		

Package

100 ug

Storage

Store at -20°C for up to 1 year. After opening, store at 2-8°C. Protect from light and moisture. Reagents must be used within 1 week of opening.


Buffer Compatibility

Buffer Components & Conditions	
Purified antibody	Yes
Antibody in ascites fluid, serum, hybridoma or tissue culture media	No
Antibody concentration	0.5-2 mg/mL
pH	6.5-8.0
Amine free buffer (e.g. MES, MOPS, HEPES, PBS)	Yes
Non-buffering salts (e.g. sodium chloride)	Yes
BSA	Yes
Sodium Azide	Yes
Chelating agents (e.g. EDTA)	Yes
Glycerol	<10%
Sugars	Yes
Gelatin	Unknown
Tris	Yes
Glycine	No
Thiomersal / Thimerosal	Yes
Merthiolate	Yes
Proclin	Yes
Borate buffer	Yes
Primary amines (e.g. lysine, peptide or ethanolamine)	No
Heavy metal ions (Zn ²⁺ , Cu ²⁺ , Mn ²⁺ , or Fe ³⁺)	No

Required Materials

Orbital Thermal Mixer, Pipettes, Sterile PBS.

Validation Data



1x10⁶ mouse bone marrow cells were surface stained with FcZero-rAb® CoraLite® Plus 488 Anti-Mouse CD11b and 0.1 ug APC-Cyanine7 Anti-Mouse Ly-6G Rabbit RecAb (AY7-98284, Clone: 242141B11) or APC-Cyanine7 Rabbit IgG Isotype Control Recombinant Antibody (AY7-98136, Clone: 240953C9). Cells were incubated with MonoZero™... Monocytes blocking Reagent (PF00020) prior to staining. Cells were not fixed.