

For Research Use Only

UltraPolymer Goat anti-Rabbit IgG (H&L)-HRP

Catalog Number: PR30011



Description

UltraPolymer Goat anti-Rabbit IgG (H&L)-HRP is a polymer (secondary antibody) formed by coupling multiple peroxidase (HRP) molecules and anti rabbit immunoglobulin (IgG) onto the polymer using the latest polymerization technology. The polymer (secondary antibody) has the characteristics of amplifying signals to increase sensitivity and reducing non-specific background. It can be used together with rabbit first antibody, and then react with Chromogen(DAB) to obtain brown color.

Package

5 mL


Storage

Store the reagent at 2-8°C. It is stable for 6 months from the date of receipt.

Application Note

This reagent is recommended for FFPE tissues and optimized for manual staining.

Validation Data



Immunohistochemistry analysis of paraffin-embedded human kidney tissue slide using PR30011 without (left) or with (right) primary antibody.

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB(1-888-478-4522)(toll free in USA), or
1(312) 455-8498(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.