

For Research Use Only

FcZero-rAb® PE-Cyanine7 Anti-Rat CD11b/c (OX42) Rabbit Recombinant Antibody

Catalog Number: PY7-FcA65606



Basic Information

Catalog Number:

PY7-FcA65606

Size:

100ug , 200 ug/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_012711

GeneID (NCBI):

25021

UNIPROT ID:

Q9J130

Full Name:

integrin alpha M

Calculated MW:

127 kDa

Purification Method:

Protein A purification

CloneNo.:

OX42

Recommended Dilutions:

FC: 0.1 ug per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

496 nm, 565 nm / 778 nm

Applications

Tested Applications:

FC

Species Specificity:

rat

Positive Controls:

FC : rat bone marrow cells,

Background Information

Integrins are cell adhesion receptors that are heterodimers composed of non-covalently associated α and β subunits (PMID: 9779984). CD11b (integrin alpha M or CR3A) and CD11c (integrin alpha X) are members of the integrin alpha chain family and can form α/β heterodimers with CD18 (integrin β 2). CD11b/CD18 is a receptor for the complement protein fragment C3bi, and is also a receptor for fibrinogen, factor X, and ICAM1 (PMID: 2971974; 15485828). CD11c/CD18 acts as a receptor for fibrinogen and is important in monocyte adhesion and chemotaxis (PMID: 1671533). The OX-42 antibody recognizes a common epitope shared between CD11b and CD11c. It binds to C3bi receptor and is known to block C3bi binding (PMID: 3512425; 8449218).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ rat bone marrow cells were surface stained with 0.1 ug PE-Cyanine7 Anti-Rat CD11b/c (OX42) Rabbit IgG RecAb (PY7-FcA65606, Clone: OX42) or FcZero-rAb™ PE-Cyanine7 Rabbit IgG Isotype Control Recombinant Antibody (PY7-FcA98136, Clone: 240953C9). Cells were incubated with MonoZero™ Monocytes blocking Reagent (PF00020) prior to staining. Cells were not fixed.



1x10⁶ rat bone marrow cells were surface stained with 0.1 ug PE-Cyanine7 Anti-Rat CD11b/c (OX42) Rabbit IgG RecAb (PY7-FcA65606, Clone: OX42) (red) or FcZero-rAb™ PE-Cyanine7 Rabbit IgG Isotype Control Recombinant Antibody (PY7-FcA98136, Clone: 240953C9) (blue). Cells were incubated with MonoZero™ Monocytes blocking Reagent (PF00020) prior to staining. Cells were not fixed.