

For Research Use Only

FcZero-rAb® PE-Cyanine7 Anti-Human Tissue Factor/CD142 Rabbit Recombinant Antibody

Catalog Number: PY7-FcA98104



Basic Information

Catalog Number:

PY7-FcA98104

Size:

100 tests, 5 ul/test

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1052

GenBank Accession Number:

NM_001993.5

GeneID (NCBI):

2152

UNIPROT ID:

P13726-1

Full Name:

coagulation factor III (thromboplastin, tissue factor)

Calculated MW:

33 kDa

Purification Method:

Protein A purification

CloneNo.:

240784E2

Recommended Dilutions:

FC: 5 ul per 10⁶ cells in a 100 µl suspension

Excitation/Emission maximum wavelengths:

496 nm, 565 nm / 778 nm

Applications

Tested Applications:

FC

Species Specificity:

human

Positive Controls:

FC : LPS treated human PBMCs,

Background Information

Normal blood coagulation is a complex process, involving a cascade of activation of different plasma proteins, ultimately resulting in the formation of a clot, called fibrin. Tissue Factor (TF), also named CD142, coagulation factor III or F3, is the primary initiator of the blood coagulation cascade and plays an essential role in hemostasis (PMID: 26877187; PMID: 33796108).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3



For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ untreated or LPS-treated human PBMCs were surface stained with FcZero-rAb® CoraLite® Plus 405 Anti-Human CD14 and 5 ul PE-Cyanine7 Anti-Human Tissue Factor/CD142 Rabbit RecAb (PY7-FcA98104, Clone: 240784E2). Cells were incubated with MonoZero™ Monocytes blocking Reagent (PF00020) prior to staining. Cells were not fixed.

1x10⁶ LPS-treated human PBMCs were surface stained with 5 ul PE-Cyanine7 Anti-Human Tissue Factor/CD142 Rabbit RecAb (PY7-FcA98104, Clone: 240784E2) (red) or FcZero-rAb® PE-Cyanine7 Rabbit IgG Isotype Control Recombinant Antibody (PY7-FcA98136, Clone: 240953C9) (blue). 1x10⁶ untreated human PBMCs were surface stained with 5 ul PE-Cyanine7 Anti-Human Tissue Factor/CD142 Rabbit RecAb (PY7-FcA98104, Clone: 240784E2) (black, dashed).

