

For Research Use Only.
Not For Use In Diagnostics.

Multi-rAb® Biotin-Goat Anti-Mouse Recombinant Secondary Antibody (H+L)



Catalog Number: **RGAM009**

Information

Catalog Number:

RGAM009

Host:

Goat

Applications:

WB, IHC, ELISA

Reactivity:

Mouse

Physical State:

Liquid

Conjugation:

Biotin

Clonality:

Multiclonal recombinant

Recommended Dilutions

1:100-1:1000 for enzyme immunohistochemistry on tissue sections;
1:5000-1:20000 for ELISA and Western blotting with chromogenic substrates.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 10 mg/mL BSA, 0.1% ProClin 300, pH 7.4.

Aliquoting is unnecessary

Purity

The antibody was purified from culture media supernatant by immunoaffinity chromatography using Protein G beads.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human ovary cancer tissue slide using 66379-1-Ig (Beta Catenin antibody) at dilution of 1:5000. Multi-rAb® Biotin-Goat Anti-Mouse Recombinant Secondary Antibody (H+L) was used at 6.5 ug/mL for signal amplification. HRP conjugated streptavidin was used for detection.