

For Research Use Only.  
Not For Use In Diagnostics.

# Multi-rAb® Biotin-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L)



Catalog Number: **RGAR009**

## Information

|                                 |                           |                                       |
|---------------------------------|---------------------------|---------------------------------------|
| Catalog Number:<br>RGAR009      | Reactivity:<br>Rabbit     | Clonality:<br>Multiclonal recombinant |
| Host:<br>Goat                   | Physical State:<br>Liquid |                                       |
| Applications:<br>WB, IHC, ELISA | Conjugation:<br>Biotin    |                                       |

## Recommended Dilutions

1:100-1:1000 for enzyme immunohistochemistry on tissue sections;  
1:5000-1:20000 for ELISA and Western blotting with chromogenic substrates.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 50% glycerol, 10 mg/mL BSA, 0.1% ProClin 300, pH 7.4.  
**Aliquoting is unnecessary**

## Purity

The antibody was purified from culture media supernatant by immunoaffinity chromatography using Protein G beads.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human appendicitis tissue slide using 80008-1-RR (Alpha smooth muscle actin specific antibody) at dilution of 1:10000. Multi-rAb® Biotin-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) was used at 4 ug/mL for signal amplification. HRP conjugated streptavidin was used for detection.