

For Research Use Only.  
Not For Use In Diagnostics.

# HRP-Goat Anti-Rabbit IgG Conformation Specific Recombinant Secondary Antibody



Catalog Number: **RGAR301**

## Information

Catalog Number:

RGAR301

Host:

Goat

Applications:

IP-WB, WB, ELISA

Reactivity:

Rabbit

Physical State:

Liquid

Conjugation:

HRP

Clonality:

Recombinant monoclonal

## Recommended Dilutions

1:3000-1:10,000 for ELISA;

1:2500-1:5000 for western blotting or IP-WB with ECL substrates.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 10 mg/mL BSA, 0.1% Proclin-300, pH 7.4.

Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from culture media supernatant by immunoaffinity chromatography using Protein G beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

