

For Research Use Only

Multi-rAb® Cyclin B1 Multi-Recombinant antibody

Catalog Number: RMX00076



Basic Information

Catalog Number: RMX00076	GenBank Accession Number: BC006510	Purification Method: N/A
Size: 100ul , Concentration: 1040 µg/ml by Nanodrop;	GeneID (NCBI): 891	Recommended Dilutions: WB: 1:5000-1:50000
Source: Rabbit	UNIPROT ID: P14635	
Isotype: IgG	Full Name: cyclin B1	
	Calculated MW: 48 kDa	
	Observed MW: 55-60 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HT-29 cells, HeLa cells
Species Specificity: human	

Background Information

Cyclin B1 is a regulatory protein involved in mitosis. The gene product complexes with p34(cdc2) to form the maturation-promoting factor (MPF). Two alternative transcripts have been found, a constitutively expressed transcript and a cell cycle-regulated transcript, that is expressed predominantly during G2/M phase of the cell cycle. The different transcripts result from the use of alternate transcription initiation sites. The antibody is specific to CCNB1. We got a 55-60 kDa band in western blotting maybe due to phosphorylation.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with RMX00076 (Cyclin B1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.