

For Research Use Only.
Not For Use In Diagnostics.

HRP-conjugated Affinipure Donkey Anti-Rat IgG(H+L)

Catalog Number:SA00001-10



www.ptglab.com

Information

Catalog Number:
SA00001-10
Host:
Donkey
Applications:
ELISA, WB, Dot blot

Reactivity:
Rat (Minimal Cross Reactivity:
Bovine, Chicken, Goat, Guinea Pig,
Syrian Hamster, Horse, Human,
Mouse, Rabbit, Sheep Serum
Proteins)
Physical State:
Liquid
Conjugation:
HRP

Recommended Dilutions

1:2000-1:10,000 for western blotting with ECL substrates
1:1000-1:20,000 for ELISA and Western blotting with chromogenic substrates

Stabilizer

3 mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.
Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data