

For Research Use Only.
Not For Use In Diagnostics.

HRP -Protein A

Catalog Number:SA00001-18



www.ptglab.com

Information

Catalog Number:

SA00001-18

Host:

Applications:

WB, Dot blot

Reactivity:

FC Fragment

Physical State:

Liquid

Conjugation:

HRP

Recommended Dilutions

1:4000-1:10000 for Western Blotting with chromogenic substrates after immunoprecipitation

Stabilizer

3 mg/ml BSA

Safety Notes

1.This product is for research use only,not for diagnostic or therapeutic use;
2. Do not use Sodium Azide as preservative.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.

Aliquoting is unnecessary for -20°C storage

Purity

Protein A was purified by chromatography from Staphylococcus aureus


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP Result of anti-EBAG9 (IP:12255-1-AP, 3ug;
Detection:12255-1-AP 1:800) with mouse heart
tissue lysate. 5A00001-18 (HRP-conjugated Protein
A from *S.aureus*) as secondary antibody.