

For Research Use Only.
Not For Use In Diagnostics.

HRP-conjugated Recombinant Rabbit Anti-Mouse IgG Kappa Light Chain



Catalog Number: **SA00001-19**

Information

Catalog Number:
SA00001-19
Host:
Rabbit
Applications:
ELISA, WB, Dot blot

Reactivity:
Mouse
Physical State:
Liquid
Conjugation:
HRP

Clonality:
Recombinant

Recommended Dilutions

1:2000-1:20000

Preservative

Proclin300

Stabilizer

0.5% BSA

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.03% Proclin300, 50% glycerol, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

Purity

Protein A purification


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Mouse serum lysates were subjected to SDS PAGE followed by western blot with SA00001-19 (mouse IgG Kappa Light Chain antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.