

For Research Use Only.  
Not For Use In Diagnostics.

# Alkaline Phosphatase-conjugated Donkey Anti-Rabbit IgG(H+L)

Catalog Number:SA00002-5



## Information

### Catalog Number:

SA00002-5

### Host:

Donkey

### Applications:

ELISA, WB

### Reactivity:

Rabbit

### Physical State:

Liquid

### Conjugation:

Alkaline phosphatase

### Clonality:

Polyclonal

## Recommended Dilutions

1:1000 – 1:3000 for ELISA and Western blotting

## Preservative

0.01% Sodium Azide

## Stabilizer

5 mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

0.01 M Tris-HCl, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.01% Sodium Azide, pH 8.0.

Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

