

For Research Use Only.  
Not For Use In Diagnostics.

# Alkaline Phosphatase–conjugated Affinipure Donkey Anti-Rat IgG(H+L)

Catalog Number:SA00002-6



## Information

Catalog Number:	SA00002-6	Reactivity:	Rat
Host:	Donkey	Physical State:	Liquid
Applications:	ELISA, WB	Conjugation:	Alkaline phosphatase

## Recommended Dilutions

1:1000 – 1:3000 for ELISA and Western blotting

## Preservative

0.01% Sodium Azide

## Stabilizer

5mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
0.01 M Tris-HCl, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.01% Sodium Azide, pH 8.0.  
Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data