

For Research Use Only.  
Not For Use In Diagnostics.

# Alkaline Phosphatase–conjugated Affinipure Goat Anti-Human IgG(H+L)



[www.ptglab.com](http://www.ptglab.com)

Catalog Number: SA00002-8

## Information

|                 |           |                 |                      |
|-----------------|-----------|-----------------|----------------------|
| Catalog Number: | SA00002-8 | Reactivity:     | Human                |
| Host:           | Goat      | Physical State: | Liquid               |
| Applications:   | ELISA, WB | Conjugation:    | Alkaline phosphatase |

## Recommended Dilutions

1:1000 – 1:3000 for ELISA and Western blotting

## Preservative

0.01% Sodium Azide

## Stabilizer

5mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
0.01 M Tris-HCl, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.01% Sodium Azide, pH 8.0.  
Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data