

For Research Use Only.
Not For Use In Diagnostics.

Fluorescein (FITC)-conjugated Affinipure Donkey Anti-Human IgG(H+L)



www.ptglab.com

Catalog Number: SA00003-10

Information

Catalog Number:	SA00003-10	Reactivity:	Human
Host:	Donkey	Physical State:	Liquid
Applications:	IF, FC	Conjugation:	FITC

Recommended Dilutions

1:20-1:100 for most applications

Preservative

0.02% Sodium Azide

Fluorophore

Fluorescein-5-isothiocyanate (FITC) Amax=492nm; Emax=520nm

Fluorophore Protein

7.7ug/mg; 2.9 moles FITC per mole IgG

Stabilizer

5mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6
Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data