

For Research Use Only.  
Not For Use In Diagnostics.

# Fluorescein (FITC)-conjugated Affinipure Goat Anti-Rabbit IgG(H+L)



[www.ptglab.com](http://www.ptglab.com)

Catalog Number: SA00003-2

## Information

Catalog Number:

SA00003-2

Host:

Goat

Applications:

IF, FC

Reactivity:

Rabbit

Physical State:

Liquid

Conjugation:

FITC

Clonality:

Polyclonal

## Recommended Dilutions

1:100-1:1000 for most applications

## Preservative

0.02% Sodium Azide

## Fluorophore

Fluorescein-5-isothiocyanate (FITC)  $\lambda_{max}=492nm$ ;  $\lambda_{em}=520nm$

## Fluorophore Protein

9.0 ug/mg; 3.4 moles FITC per mole IgG

## Stabilizer

5 mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

Storage:

Store at  $-20^{\circ}C$ . Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for  $-20^{\circ}C$  storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data