

For Research Use Only.
Not For Use In Diagnostics.

Biotin-conjugated Affinipure Donkey Anti-Goat IgG(H+L)

Catalog Number:SA00004-3



www.ptglab.com

Information

Catalog Number:

SA00004-3

Host:

Donkey

Applications:

ELISA, WB, IHC/ICC

Reactivity:

Goat

Physical State:

Liquid

Conjugation:

Biotin

Recommended Dilutions

1:100 – 1:1000 for enzyme immunohistochemistry on tissue sections

1:4000-1:80000 for ELISA and Western blotting with chromogenic substrates

Preservative

0.01% Sodium Azide

Stabilizer

5mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.01% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data