

For Research Use Only.
Not For Use In Diagnostics.

AMCA-conjugated Goat Anti-Mouse IgG(H+L)

Catalog Number: SA00010-1



www.ptglab.com

Information

Catalog Number:

SA00010-1

Host:

Goat

Applications:

IF, FC

Reactivity:

Mouse

Physical State:

Liquid

Conjugation:

AMCA

Clonality:

Polyclonal

Recommended Dilutions

1:20-1:100 for most applications

Preservative

0.02% Sodium Azide

Fluorophore

7-Amino-4-methylcoumarin-3-acetic acid $\lambda_{max}=350nm$, $\lambda_{exc}=450nm$

Stabilizer

5 mg/ml BSA

Absorbance Ratio

A350/A280=0.76

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:

Store at $-20^{\circ}C$. Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for $-20^{\circ}C$ storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Methanol) fixed MDCK cells using Beta Actin antibody (66009-1-Ig, Clone: 2D4H5) at dilution of 1:400 and AMCA-conjugated Goat Anti-Mouse IgG(H+L)(SA00010-1).