

For Research Use Only.  
Not For Use In Diagnostics.

# AMCA-conjugated Donkey Anti-Goat IgG(H+L)

Catalog Number: SA00010-3



[www.ptglab.com](http://www.ptglab.com)

## Information

### Catalog Number:

SA00010-3

### Host:

Donkey

### Applications:

IF, FC

### Reactivity:

Goat

### Physical State:

Liquid

### Conjugation:

AMCA

### Clonality:

Polyclonal

## Recommended Dilutions

1:20-1:100 for most applications

## Preservative

0.02% Sodium Azide

## Fluorophore

7-Amino-4-methylcoumarin-3-acetic acid  $\lambda_{max}=350nm$ ,  $\lambda_{em}=450nm$

## Stabilizer

5 mg/ml BSA

## Absorbance Ratio

A350/A280=0.71

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

### Storage:

Store at  $-20^{\circ}C$ . Stable for one year after shipment.

### Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for  $-20^{\circ}C$  storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

