

For Research Use Only.
Not For Use In Diagnostics.

CoraLite594 – conjugated Donkey Anti-Mouse IgG(H+L)

Catalog Number:SA00013-7



www.ptglab.com

Information

Catalog Number:

SA00013-7

Host:

Donkey

Applications:

IF, FC

Reactivity:

Mouse (Minimal Cross Reactivity:

Bovine, Chicken, Goat, Guinea Pig,

Syrian Hamster, Horse, Human,

Rabbit, Sheep Serum Proteins)

Physical State:

Liquid

Conjugation:

CoraLite®594

Clonality:

Polyclonal

Recommended Dilutions

1:100 – 1:500

Preservative

0.02% Sodium Azide

Fluorophore

CL594 Fluorescent Dye
A_{max}=593 nm, E_{max}=614 nm

Stabilizer

5 mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Methanol) fixed MDCK cells using Beta Actin antibody (66009-1-Ig, Clone: 2D4H5) at dilution of 1:400 and CoraLite594 - conjugated Donkey Anti-Mouse IgG(H+L)(SA00013-7).