

For Research Use Only.
Not For Use In Diagnostics.

CoraLite488-conjugated Recombinant Rabbit Anti-Mouse IgG Kappa Light Chain



Catalog Number:SA00014-3

Information

Catalog Number:

SA00014-3

Host:

Rabbit

Applications:

IF

Reactivity:

Mouse

Physical State:

Liquid

Conjugation:

CoraLite®488

Clonality:

Recombinant

Recommended Dilutions

1:100-1:500

Preservative

Proclin300

Fluorophore

CL488 Fluorescent Dye, Amax=488 nm; Emax=515 nm

Stabilizer

0.5% BSA

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

Purity

Protein A purification

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data