

For Research Use Only.  
Not For Use In Diagnostics.

# CoraLite594-conjugated Mouse Anti-Heavy Chain of Rabbit IgG

Catalog Number:SA00014-5



## Information

### Catalog Number:

SA00014-5

### Host:

Mouse

### Applications:

IF, FC

### Reactivity:

Rabbit

### Physical State:

Liquid

### Conjugation:

CoraLite®594

### Clonality:

Monoclonal

## Recommended Dilutions

1:100-1:500

## Preservative

Proclin300

## Fluorophore

CL594 Fluorescent Dye. Amax=593 nm, Emax=614 nm

## Stabilizer

0.5% BSA

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

## Purity

Protein G purification

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB  
(1-888-478-4522)  
(toll free in USA),  
or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

