

For Research Use Only

anti-Mouse IgG3 Magnetic VHH Agarose for Immunoprecipitation



www.ptglab.com

Catalog Number: mIG3ma

Basic Information

Catalog Number:
mIG3ma

Applications:
IP, Co-IP

Conjugate:
Magnetic Agarose beads; ~40 um (cross-linked 6% magnetic agarose beads)

Host:
Alpaca

Type:
Nanobody

Class:
Recombinant - Animal free production

Description

anti-Mouse IgG3 IP Beads is an affinity resin for IP of Mouse IgG3. It consists of rabbit specific VHHs (Nanobodies) coupled to Magnetic Agarose beads.

Binding capacity

/

Elution buffer

SDS Sample Buffer

Wash buffer compatibility

/

Affinity (K_D)

/

Storage

Storage:
+4°C / do not freeze!

Storage Buffer:
20% Ethanol

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.