

Catalog Number: mta

2 Publications

Basic Information

Catalog Number:
mta

Applications:
IP, CoIP

Conjugate:
Agarose beads; bead size: ~ 90 µm (cross-linked 4 % agarose beads)

Host:
Alpaca

Type:
Nanobody

Class:
Recombinant

Description

ChromoTek's MK2-Trap Agarose are affinity beads for IP of MK2-fusion proteins. It comprises a MK2 Nanobody/ VHH conjugated to agarose beads.

Specificity/Target

MK2 (MAP kinase-activated protein kinase 2, MAPKAP-K2)

Elution buffer

SDS sample buffer
0.2 M glycine pH 2.5

Storage

Storage:
Upon receipt store at +4°C. Do not freeze!

Storage Buffer:
20% ethanol

Selected Validation Data

