

For Research Use Only

anti-mNeonGreen recombinant VHH, unconjugated



www.ptglab.com

Catalog Number: nt

1 Publications

Basic Information

Catalog Number:

nt

Applications:

Suitable for conjugation of dyes, biotin, beads, surfaces etc. via NHS ester reaction

Host:

Alpaca

Conjugate:

Unconjugated

Type :

Nanobody

Class:

Recombinant

RRID:

AB_2827558

Molecular Weight:

14.2 kDa

Description

Alpaca anti-mNeonGreen VHH, purified recombinant binding protein

Specificity/Target

mNeonGreen of Branchiostoma lanceolatum

For the complete list, please click here: [Fluorescent protein specificity table](#)

Affinity (K_D)

Dissociation constant K_D of 2 nM

Storage

Storage:

Shipped at ambient temperature. Upon receipt store at 4°C; stable for one year. Do not freeze.

Storage Buffer:

25mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, preservative: 0.09 % sodium azide

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under
Proteintech Group brand and is not available
to purchase from any other manufacturer.

Selected Validation Data



Structure of mNeonGreen.