

For Research Use Only

iST mNeonGreen-Trap IP-MS Sample Preparation Kit



www.ptglab.com

Catalog Number: ntak-iST

1 Publications

Basic Information

Catalog Number:
ntak-iST

Applications:
IP, MS

Conjugate:
Agarose beads; bead size: ~ 90 µm (cross-linked 4 % agarose beads)

Host:
Alpaca

Type:
Nanobody

Class:
Recombinant

Description

Kit for IP/ Co-IP of mNeonGreen tagged proteins and sample preparation for mass spectrometry (MS)

Binding capacity

30 µg of recombinant mNeonGreen per 25 µL bead slurry

Specificity/Target

mNeonGreen

Affinity (K_D)

Dissociation constant K_D of 2 nM

Storage

Storage:

Kit is shipped at ambient temperature. Upon receipt store mNeonGreen-Trap at 4°C and Digest at -20°C


Storage Buffer:
20% EtOH

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Schematic workflow of iST Nano-Trap® Kit

Save up to 40 hours on your mass spec experiments with the iST Nano-Trap® Kits