

For Research Use Only

Flow Cytometry Staining Buffer (1X)

Catalog Number: PF00012



www.ptglab.com

Description

Ready-to-use Flow Staining Buffer (1X) is recommended immunofluorescent staining protocols using viable or fixed cells in suspension for flow cytometry. Flow Staining Buffer is a buffered saline solution containing Bovine Calf Serum as a protein source to minimize nonspecific binding of antibodies and sodium azide (0.09%) as a preservative. This buffer can be used for diluting antibody, cell suspension, and as a wash steps for cell surface staining and flow cytometric analysis.

Package

500 mL

Storage

4°C

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB(1-888-478-4522)(toll free in USA), or
1(312) 455-8498(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.