

# Neural Cell Markers Essentials Antibody Kit

Catalog Number: PK30026

## Description

The Neural Cell Markers Essentials Antibody Kit provides a cost-effective tool for identifying and studying neuronal and glial cells. Perfect for researchers starting a new project, screening multiple prospective targets or those who simply require less volume.

## Product Information

The Neural Cell Markers Essentials Antibody Kit contains antibodies against 6 key proteins for identifying neuronal and glial cells.

Antigen	Catalog No.	Host, clonality	Tested Reactivity	Applications	Cell Type Marker	Volume
MAP2	<a href="#">17490-1-AP</a>	Rabbit polyclonal	H, M, Rt, Mk, G	WB, IHC, IF/IC C, IF-P, FC, IP, ELISA	Neurons (dendrites)	20 uL
NeuN	<a href="#">66836-1-Ig</a>	Mouse monoclonal	H, M, Rt, G	IHC, IF-P, FC, ELISA	Neurons (cell body)	20 uL
TUBB3	<a href="#">66375-1-Ig</a>	Mouse monoclonal	H, M, Rt, Pg, Rb, Ch	WB, IHC, IF/IC C, IF-P, FC, ELISA	Neurons (axon)	20 uL
GFAP	<a href="#">81063-1-RR</a>	Rabbit Monoclonal	H, M, Rt, Pg	WB, IHC, IF-P, ELISA	Astrocytes	20 uL
Iba1	<a href="#">81728-1-RR</a>	Rabbit Monoclonal	H, M, Rt	WB, IHC, IF-P, FC, ELISA	Microglia	20 uL
CNPase	<a href="#">66729-1-Ig</a>	Mouse monoclonal	H, M, Rt, Pg	WB, IF-P, FC, ELISA	Oligodendrocytes	20 uL

Also see our 'Neural Cell Markers Expanded Antibody Kit' on the following page

<https://www.ptglab.com/products/Neural-Cell-Markers-Expanded-Antibody-Kit-PK30027.htm>

## Package

6 × 20 uL

## Storage

Store at -20°C. Stable for one year from the date of receipt.


## Background Information

The central nervous system (CNS) is a complex network consisting of a diverse subset of cells including neurons, astrocytes, oligodendrocytes, microglia, and neural stem cells. Identifying cell populations and protein localization can facilitate studies exploring the biological processes that occur during neurodevelopment or in neurological disorders, such as cell proliferation, migration, differentiation, morphology, and synapse formation.


## Standard Protocols

Click [here](#) to view our standard protocols for various applications including WB, IP, IHC, IF, FC, and ELISA.


## Validation Data




Immunofluorescent analysis of (4% PFA) fixed rat brain tissue using MAP2 antibody (17490-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).




Immunofluorescent analysis of (4% PFA) fixed rat brain tissue using 17490-1-AP (MAP2 antibody) at dilution of 1:100 and CoraLite594-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). The section was co-stained with 60190-1-Ig (GFAP antibody, green).




Immunofluorescent analysis of (4% PFA) fixed mouse brain tissue using MAP2 antibody (17490-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).




Immunofluorescence analysis of human pluripotent stem cell-derived astrocytes with S100β (15146-1-AP) at 1/200 (Magenta) and neurons with TUJ1 (66375-1-Ig) at 1:500 (Green). The sample was fixed with 4% Paraformaldehyde and permeabilized with 0.3% Triton X-100. Alexa Fluor 488-conjugated goat anti-mouse IgG (1/500) and Alexa Fluor 594-conjugated goat anti-rabbit IgG (1/500) were used as the secondary antibodies. Nuclei were




Retinal organoids (day 30) generated from human induced pluripotent stem cells (iPSCs) and fixed with 4% PFA. Stained for Tubulin beta 3/TUJ1 using 66375-1-Ig at 1:500 dilution (green) and PAX6 (12323-1-AP) at 1:500. Nuclear stain DAPI (blue). Scale bar = 50 μm. Data generated by Alessandro Bellapianta at Johannes Kepler Universitat, Austria.




Immunofluorescent analysis of (4% PFA) fixed brain tissue using 66375-1-Ig (TUBB3-specific antibody), at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). The section was co-stained with 26975-1-AP (NeuN antibody, red).




Immunofluorescent analysis of (4% PFA) fixed rat cerebellum tissue using 66836-1-Ig (NeuN antibody, red), at dilution of 1:200 and CoraLite®594-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). The section was co-stained with 14479-1-AP (Calbindin-D28k antibody, green).




Immunofluorescent analysis of (4% PFA) fixed rat brain tissue using GFAP antibody (81063-1-RR, Clone: 4C6) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).




Immunofluorescent analysis of (4% PFA) fixed mouse brain tissue using GFAP antibody (81063-1-RR, Clone: 4C6) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunofluorescent analysis of (4% PFA) fixed rat brain tissue using IBA1 antibody (81728-1-RR, Clone: 4D5) at dilution of 1:500 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunofluorescent analysis of (4% PFA) fixed mouse brain tissue using 66729-1-Ig (CNPase antibody), at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 17490-1-AP (MAP2 antibody) at dilution of 1:5000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).

**For technical support and original validation data for this product please contact**

T: 1 (888) 4PTGLAB(1-888-478-4522)(toll free in USA), or 1(312) 455-8498(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**