

For Research Use Only

anti-p53 (C-term) recombinant VHH, unconjugated



www.ptglab.com

Catalog Number: **pt2**

Basic Information

Catalog Number:
pt2

Applications:
Suitable for conjugation of dyes, biotin, beads, surfaces etc. via NHS ester reaction

Host:
Alpaca

Conjugate:
Unconjugated

Type :
Nanobody

Class:
Recombinant

RRID:
AB_2631381

Molecular Weight:
14.8 kDa

Description

Alpaca anti-p53 VHH, purified recombinant binding protein

Specificity/Target

p53
Species-Reactivity: Tested on human, mouse
Isoform: Alpha
Epitope: 302-393 aa

Storage

Storage:
Shipped at ambient temperature. Upon receipt store at 4°C; stable for one year. Do not freeze.

Storage Buffer:
25mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, preservative: 0.09 % sodium azide

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under
Proteintech Group brand and is not available
to purchase from any other manufacturer.

Selected Validation Data



Structure of VHH.