

Catalog Number: pta

9 Publications

Basic Information

Catalog Number:

pta

Applications:

IP, CoIP

Conjugate:

Agarose beads; bead size: ~ 90 µm (cross-linked 4 % agarose beads)

Host:

Alpaca

Type:

Nanobody

Class:

Recombinant

Description

p53-N-term-Trap Agarose is an affinity reagent for IP of p53. It consists of a p53 Nanobody/ VHH coupled to agarose beads, and which binds to the N-terminus of p53.

Specificity/Target

p53 (isoforms: alpha, beta, and gamma); epitope 1-81 aa

Elution buffer

SDS sample buffer
0.2 M glycine pH 2.5

Storage


Storage:

Upon receipt store at +4°C. Do not freeze!

Storage Buffer:

20% ethanol

Selected Validation Data



p53-N-term-Trap for immunoprecipitation of p53 isoform alpha, beta and gamma. I: Input, FT: Flow-through, B: Bound.