

For Research Use Only

Nano-Secondary® anti-human IgG/anti-rabbit IgG, recombinant VHH, for 2x Cys conjugation [CTK0101]

Catalog Number: **srbGCys2-1**

2 Publications



www.ptglab.com

Basic Information

Catalog Number:

srbGCys2-1

Applications:

IF, WB, FC

Host:

Alpaca

Conjugate:

N- and C-terminal cysteine conjugation with thiol-reactive reagents

Type:

Nano-Secondary®, secondary Nanobody (VHH)

RRID:

AB_2864261

Description

Monovalent, recombinant secondary single domain antibody to human IgG and rabbit IgG: alpaca monoclonal Nanobody, Fc-specific, for cysteine conjugation with thiol-reactive reagents, e.g. maleimides

Specificity/Target

Fc-fragment of human IgG and rabbit IgG (co-reactivity)

Affinity (K_D)

K_D = 0.2 nM

Storage

Storage:

Shipped on dry ice. Store at -20°C/-4°F. Avoid freeze-thaw cycles.

Storage Buffer:

10 mM HEPES pH 7.0, 500 mM NaCl, 1 mM TCEP
Preservative: 0.09 % Sodium azide

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.