

For Research Use Only

anti-GST recombinant VHH, unconjugated



www.ptglab.com

Catalog Number: st

4 Publications

Basic Information

Catalog Number:

st

Applications:

Suitable for conjugation of dyes, biotin, beads, surfaces etc. via NHS ester reaction

Host:

Alpaca

Conjugate:

Unconjugated

Type :

Nanobody

Class:

Recombinant

RRID:

AB_2631368

Molecular Weight:

14.7 kDa

Description

Alpaca anti-GST VHH, purified recombinant binding protein

Specificity/Target

Glutathione S-Transferase (GST) from Schistosoma japonicum

Affinity (K_D)

Dissociation constant K_D of 1 nM

Storage

Storage:

Shipped at room temperature. Upon receipt store at 4°C; stable for one year. Do not freeze.

Storage Buffer:


25mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, preservative: 0.09 % sodium azide

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Structure of GST.