

TurboGFP recombinant antibody, VHH-rabbit IgG Fc fusion [CTK0204]



Product code: tbrb

Properties

Description	Rabbit monoclonal antibody to TurboGFP. Recombinant fusion of anti-TurboGFP Nanobody (VHH) [CTK0204] and rabbit IgG Fc domain.
Target / Specificity	TurboGFP
Species Specificity	Recombinant protein from <i>Pontellina plumata</i>
GenBank Accession Number	AY268072
GeneID (NCBI)	N/a
Product Type	Primary antibody; chimeric fusion of Nanobody (VHH) and rabbit IgG Fc domain
Format	Recombinant monoclonal rabbit antibody
Source	Rabbit, recombinantly produced
Isotype / Subclass	Rabbit IgG (CH2 and CH3 domain of rabbit IgG)
Clonality	Monoclonal
Clone	CTK0204
RRID	AB_2889371
Immunogen	Purified recombinant TurboGFP [<i>Pontellina plumata</i>]
Conjugate	Unconjugated
Purification Method	Recombinantly expressed and purified by IMAC. <i>Note: Recombinant VHH-Fc fusion antibody carries a His-tag.</i>
Form	Buffered aqueous solution
Concentration	1000 µg/mL
Calculated MW	70.0 kDa
Tested Applications	IF
Positive Control	IF: HeLa cells transfected with Actin-Chromobody®-TurboGFP
Cited Applications	
Recommended Dilutions	IF/ICC: 1:400 - 1:2,000
Storage Buffer	PBS, 10% glycerol Preservative: 0.02 % sodium azide, safety datasheet (SDS): sodium azide
Storage Conditions	Aliquot upon receipt and store at -20°C/-4°F. Avoid freeze-thaw cycles.
Stability	Stable for 1 year at -20°C/-4°F after shipment.
Shipment	Shipped on ice

TurboGFP recombinant antibody, VHH-rabbit IgG Fc fusion [CTK0204]

Product code: tbrb


Size	20 µL; 100 µL
Synonyms	VHH, Nanobody, alpaca single domain antibody, binding domain of single domain antibody, Nano-antibody

Selected validation data

Immunofluorescence

Spin down material before use!

Dilute the antibody 1:400 - 1:2,000 in blocking buffer (e.g. PBS supplemented with 4% BSA) and incubate for 1 h at room temperature.



HeLa cells transiently expressing Actin-Chromobody®-TurboGFP were immunostained with TurboGFP recombinant antibody, VHH-rabbit IgG Fc fusion [CTK0204] (tbfrb-20, 1:1,000) and Nano-Secondary® alpaca anti-human IgG/anti-rabbit IgG, recombinant VHH, Alexa Fluor® 647 [CTK0101, CTK0102] (srbAF647-1, 1:1,000). Scale bar, 10 µM.

TurboGFP recombinant antibody, VHH-rabbit IgG Fc fusion [CTK0204]



Product code: tbrb

Product overview and related products

TurboGFP toolbox	Product code
TurboGFP-Trap Agarose	tbta-10; -20; -100
TurboGFP-Trap Agarose Kit	tbtak-20
TurboGFP-Trap Magnetic Agarose	tbtma-10; -20; -100
TurboGFP-Trap Magnetic Agarose Kit	tbtmak-20
IST TurboGFPn-Trap Kit (for IP/MS)	tbtak-IST-8
Binding Control Agarose	bab-20
Binding Control Magnetic Agarose	bmab-20
Spin columns	sct-10; sct-20; sct-50
TurboGFP VHH, recombinant binding protein	tbt-250
TurboGFP recombinant antibody, VHH-mouse IgG1 Fc fusion [CTK0204]	tbfms-20; -100
TurboGFP recombinant antibody, VHH-rabbit IgG Fc fusion [CTK0204]	tbfrb-20; -100
Nano-Secondary [®] alpaca anti-mouse IgG1, recombinant VHH, Alexa Fluor [®] 488 [CTK0103, CTK0104]	sms1AF488-1-10; -100
Nano-Secondary [®] alpaca anti-mouse IgG1, recombinant VHH, Alexa Fluor [®] 568 [CTK0103, CTK0104]	sms1AF568-1-10; -100
Nano-Secondary [®] alpaca anti-mouse IgG1, recombinant VHH, Alexa Fluor [®] 647 [CTK0103, CTK0104]	sms1AF647-1-10; -100
Nano-Secondary [®] alpaca anti-human IgG/anti-rabbit IgG, recombinant VHH, Alexa Fluor [®] 488 [CTK0101, CTK0102]	srbAF488-1-10; -100
Nano-Secondary [®] alpaca anti-human IgG/anti-rabbit IgG, recombinant VHH, Alexa Fluor [®] 568 [CTK0101, CTK0102]	srbAF568-1-10; -100
Nano-Secondary [®] alpaca anti-human IgG/anti-rabbit IgG, recombinant VHH, Alexa Fluor [®] 647 [CTK0101, CTK0102]	srbAF647-1-10; -100

For product details, information, and ordering visit www.chromotek.com and www.ptglab.com.

TurboGFP recombinant antibody, VHH-rabbit IgG Fc fusion [CTK0204]



Product code: tbrb

Contact

support@chromotek.com

ChromoTek GmbH
Am Klopferspitz 19
82152 Planegg-Martinsried
Germany
phone: +49 89 124 148 80
fax: +49 89 124 148 811

ChromoTek Inc.
62-64 Enter Lane
Islandia, NY 11749
USA
phone: 631 501 1058
fax: 631 501 1060

Disclaimer

Only for research applications, not for diagnostic or therapeutic use!

The VHH-Fc fusion recombinant antibodies have been developed in collaboration with Absolute Antibody Ltd.

ChromoTek and GFP-Trap, RFP-Trap, Myc-Trap, Spot-Trap, Spot-Tag, Spot-Label, Spot-Cap, Nano-Secondary, F2H Kit, and Chromobody are registered trademarks of ChromoTek GmbH, part of Proteintech Group. Nano-CaptureLigand and V5-Trap are trademarks of ChromoTek GmbH, part of Proteintech Group. Nanobody is a registered trademark of Ablynx, a Sanofi company. Alexa Fluor is a registered trademark of Life Technologies Corporation, a part of Thermo Fisher Scientific Inc. Dynabeads is a trademark of Life Technologies AS, a part of Thermo Fisher Scientific Inc. SNAP-tag is a registered trademark and CLIP-tag is a trademark of New England Biolabs, Inc. Octet is a registered trademark of FortéBio, a Sartorius brand. Other suppliers' products may be trademarks or registered trademarks of the corresponding supplier each. Statements on other suppliers' products are given according to our best knowledge.